

Summer Sale on Nanodiamond Products

Save Up to 50%

(July 1st - August 15th 2017)

DETONATION NANODIAMOND



*Fully deagglomerated, 4-5 nm primary particles of detonation nanodiamond with negative surface charge suspended in deionized water (DI water) at 1% (w/v) (10 g/L)

Product	Catalog #	Regular price	Sale price
Monodispersed 5 nm DND, 500 mL	ND5nmNH2O500ml	\$400	\$300 (25% off)
Monodispersed 5 nm DND, 1000 mL	ND5nmNH2O1000ml	\$750	\$450 (40% off)
Average particle size 10-15 nm, 500 mL	ND10nmNH2O500ml	\$375	\$280 (25% off)
Average particle size 10-15 nm, 100 mL	ND10nmNH2O1000ml	\$600	\$360 (40% off)

**Maximum volume of purchase is limited to 10 L per customer*

FLUORESCENT NANODIAMOND



Fluorescent, high-pressure high-temperature (HPHT) synthetic nanodiamonds containing Nitrogen-Vacancy (NV) centers suspended in deionized water (DI water) at 0.1% (w/v) (1 g/L). Amphoteric surface chemistries with negative charge.

High Brightness (~3 ppm NV centers)

Product	Catalog #	Regular price	Sale price
140 nm High Brightness (1 mg/ml, 10 mL)	NDNV140nmHi10ml	\$270	\$190 (30% off)
100 nm High Brightness (1 mg/ml, 10 mL)	NDNV100nmHi10ml	\$350	\$315 (10% off)
90 nm High Brightness (1 mg/ml, 10 mL)	NDNV90nmHi10ml	\$400	\$200 (50% off)
90 nm High Brightness (1 mg/ml, 100 mL)	NDNV90nmHi100ml	\$3500	\$1750 (50% off)

Summer Sale on Nanodiamond Products

Save Up to 50%

(July 1st - August 15th 2017)

FLUORESCENT NANODIAMOND



Medium Brightness (~1.5 ppm NV centers or below)

Product	Catalog #	Regular price	Sale price
140 nm Medium Brightness(1 mg/ml, 10 mL)	NDNV140nmMd10ml	\$170	\$85 (50% off)
140 nm Medium Brightness(1 mg/ml, 100 mL)	NDNV140nmMd100ml	\$1500	\$950 (37% off)
100 nm Medium Brightness(1 mg/ml, 10 mL)	NDNV100nmMd10ml	\$240	\$120 (50% off)
100 nm Medium Brightness(1 mg/ml, 100 mL)	NDNV100nmMd100ml	\$2000	\$1000 (50% off)
60 nm Medium Brightness(1 mg/ml, 10 mL)	NDNV60nmMd10ml	\$260	\$130 (50% off)
50 nm Medium Brightness(1 mg/ml, 10 mL)	NDNV50nmMd10ml	\$260	\$130 (50% off)
40 nm Medium Brightness(1 mg/ml, 10 mL)	NDNV40nmMd10ml	\$280	\$140 (50% off)

Single Photon Emitter ND-NV (~1-4 NV per particle on average)

Product	Catalog #	Regular price	Sale price
40 nm Low Brightness(1 mg/ml, 10 mL)	NDNV40nmLw	\$150	\$110 (27% off)